

PREPARATION OF SAMPLES FOR MYCOPLASMA TESTING

1. Culture your cells for 2 weeks without any antibiotics. Split as many times during this 2 weeks as necessary just keep one flask of cells **ALWAYS** in the antibiotic free media.
2. Collect samples when cells are 80-90% confluent, **3 or more days after last split** (e.g., weekend culture) to ensure that the titer of the mycoplasmas in the medium is at a high level.
3. Collect 100 μ l of media. Boil it for 5 min at 95°C. Store at 4°C * (sample is stable for approximately 1 week).
4. Wash cells with PBS.
5. Remove adherent cells by cell scraper (NO TRYPSIN) into PBS.
6. Collect an aliquote for cell count and centrifuge cells at 4°C 1200 rpm, 5 minutes.
7. Resuspend cells in PBS, place approximately $0.5-1.0 \times 10^6$ cells into 1.5 mL Eppendorf.
8. Centrifuge at 4°C, 5000 rpm for 1 minute.
9. Remove supernatant.
10. Resuspend cell pellet with 1 mL PBS and centrifuge as per step 8.
11. Remove supernatant and store cell pellet at 4°C.*

*** The media sample as well as the cell pellet can be also stored at -20°C**